

Setting your VibraLITE 2 or VibraLITE 3 re Daylight Saving Time

When the CALENDAR CLOCK operating mode shows on the display:
Press & Hold the LAP/RESET button until the second(s) digits begin to BLINK.

Press the MODE button.
The hour digits will begin to BLINK.

Press the START/STOP button to rotate the hour until you have set the time ahead or back one hour to the correct hour (Spring forward, Fall back).

(Remember when using the 12 HR CLOCK TIME to have the PM indicator showing on the display).

Push the LAP/RESET Button and the new time is set.

If you need to change T2 (Time 2):

Press the MODE button to rotate to T2 mode.

Press & Hold the LAP/RESET button until the hour digits begin to BLINK.

Press the START/STOP button to rotate the hour until you have set the time ahead or back one hour to the correct hour (Spring forward, Fall back).

(Remember when using the 12 HR CLOCK TIME to have the PM indicator showing on the display).

Push the LAP/RESET Button and the new time is set.

Note: Your Alarm settings will be activated by the time set in the Calendar/Clock Mode.

Back to Home Page: www.globalassistive.com